

all-in-ONE™

Surgical Instrument Cleaner AND Conditioners

Highly concentrated for fast and effective cleaning.

Image: all-in-ONE™ Surgical Instrument Cleaner AND Conditioners replacing multiple products, lowering costs, and improving the performance of Surgical Instruments. The all-in-ONE™ Cleaners AND Conditioners deliver a unique formulation of detergents and multi-tiered "high-level" enzymes (Protease, Amylase, Lipase, and Carbohydrase) within a Neutral pH surfactant chemical complex. The all-in-ONE™ rapidly breaks down and removes all forms of proteinaceous bioburden, cleans and conditions the surface, enhances the "passive layer" protection of surgical stainless steel, lubricates moving parts, and then rinses clean. ONE Gallon of the all-in-ONE™ highly concentrated Surgical Instrument Cleaner AND Conditioners will deliver 512 gallons (1938 liters) of Surgical Instrument Cleaning AND Conditioning power. The Dilution rate is .25 to .5 ounces per gallon (1.92 ml per liter to 3.85 ml per liter) depending on level of bioburden, encrustation and STAINS.

The Dilution rate is .25 to .5 ounces per gallon (1.92 ml per liter to 3.85 ml per liter) depending on level of bioburden, encrustation and STAINS.

The easy FOAM-it™ lays down a thick blanket of long lasting foam to hydrate and remove proteinaceous debris, hard water deposits, and encrusted STAINS, for fast acting enzyme detergent cleaning action. When water is added the 'FOAM' dissipates providing 'free-rinsing' purity, cleaning and conditioning.

all-in-ONE™ Surgical Instrument Cleaner AND Conditioners

ONE product for

Cleaning and Conditioning Surgical Instruments.

ONE application for thoroughly

Cleaning and Conditioning Surgical Instruments.

Apply it ONCE...

It's all you really need for Cleaning AND Conditioning Surgical Instruments and Scopes.



easy FOAM-it™
Surgical Instrument
Cleaner AND Conditioners

Replace Cleaning AND Conditioning products for: removing all forms of proteinaceous bioburden, pre-cleaning and soaking to remove STAINS, manual cleaning, Ultrasonic cleaning, and automated Surgical Instrument Washing.



all-in-ONE™
Surgical Instrument
Cleaner AND Conditioners

The all-in-ONE™ Surgical Instrument Cleaners deliver a unique formulation of detergents and multi-tiered "high-level" enzymes (Protease, Amylase, Lipase, Carbohydrase) within a Neutral pH surfactant chemical complex. The all-in-ONE™ cleaners will: break down bioburden and clean the surface, remove STAINS, enhance the "passive layer" protection of surgical stainless steel, and lubricate moving parts. The highly concentrated all-in-ONE™ cleaners are aggressive against your most difficult cleaning challenge and guaranteed to be safe for all types of metal, rigid and flexible scopes, lenses, microsurgical instrumentation and reprocessing equipment.

easy FOAM-it™ Surgical Instrument Cleaning AND Conditioning 12 applicator bottles per case

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 4 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 5 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 15 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 30 gallons

all-in-ONE™ Surgical Instrument Cleaning AND Conditioning 55 gallons

Your cleaner surgical instruments will look newer and perform longer.



Product Support: Support@ergo-Logistics.com Website: www.ergo-Logistics.com
1-509-747-5027 email mail: Support@ergo-Logistics.com