

Enzymatic Surgical Instrument Detergent UPS-1 **cut costs.**

The Surgical Instrument Detergent Enzymatic Lubricants clean **residue free.**

Four Enzyme Surgical Instrument Cleaning Detergent Lubricants **clean faster.**

The Enzymatic Surgical Instrument Cleaning Detergent Lubricant product delivers:

Enzyme Cleaners (protease, amylase, carbohydrase, and lipase)

Detergent Cleaners (surfactant cleaners)

Lubricants (water soluble and free rinsing)

Surface Conditioners (strengthen the chromium passive layer of stainless steel)

Designed for use as:

manual surgical instrument and scope cleaning solution,

ultrasonic surgical instrument and scope cleaning solution,

automated surgical instrument washer and washer disinfectant concentrate.

The Enzymatic Neutral pH (7) Surgical Instrument Detergent Lubricants clean fast and residue free. The all-in-ONE Surgical Instrument Cleaners deliver the four enzymes that are necessary to break down all forms of proteinaceous bioburden and the surfactant detergents to clean the surface. The all-in-ONE Enzymatic is highly concentrated for faster cleaning and lower cleaning costs. The all-in-ONE Enzymatic delivers an average concentration that is 400% higher than the other enzyme detergents. One gallon can deliver 512 gallons (1938 liters) of enzymatic cleaning and conditioning power. The dilution rate is .25 to .5 ounces per gallon (1.9 ml liters to 3.8 ml per liter) depending on the level of bioburden, encrustation, and stains.



The all-in-ONE Enzyme Detergent Cleaner Lubricants are packaged as:

UPS-32 Enzymatic easy-FOAM-it **Pre wash** Surgical Instrument Detergent Lubricant Cleaners, 12 ea. per case,

UPS-1 Enzymatic Surgical Instrument Detergent Lubricant Cleaners, 4 gallons and 1 pump per case,

UPS-5 Enzymatic Surgical Instrument Detergent Lubricant Cleaners, 5 gallon jug,

UPS-15 Enzymatic Surgical Instrument Detergent Lubricant Cleaners, 15 gallon jug,

UPS-30 Enzymatic Surgical Instrument Detergent Lubricant Cleaners, 30 gallon jug,

UPS-55 Enzymatic Surgical Instrument Detergent Lubricant Cleaners, 55 gallon jug,



Customer Service and Purchasing: 1-509-981-1553 email: Info@yourCEBA.com

Product Technical Support: John Temple, 1-509-747-5027, email: Temple@yourCEBA.com

Websites: www.yourCEBA.com