

The UPS-1 Enzymatic Surgical Instrument Detergents **cut costs**.
Four Enzyme Surgical Instrument Cleaning Detergent Lubricants **clean faster**.
Surgical Instrument Cleaning Detergent Enzyme Lubricants clean **residue free**.

- Enzyme Cleaners (protease, amylase, carbohydrase, and lipase)
- Detergent Cleaners (surfactant cleaners)
- Lubricants (water soluble and free rinsing)
- Surface Conditioners (strengthen the chromium passive layer of stainless steel)

The Neutral pH (7) Surgical Instrument Enzymatic Detergent Lubricants clean fast and residue free.
The all-in-ONE Surgical Instrument Cleaners deliver the four enzymes that are necessary to break down all forms of proteinaceous bioburden and the surfactant detergents to clean the surface.
The all-in-ONE is highly concentrated for faster cleaning and lower cleaning costs.

The cost cutting all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners delivers the:

- highest concentration available, significantly decreasing the amount of product being used,
- reduced weight (per effective volume) reducing shipping and handling costs,
- combined enzyme and detergent cleaning treatments speeding turnaround (decontamination) times,
- lowered inventory and acquisition costs by replacing multiple products and reducing purchases,
- improved performance of surgical instruments and scopes extended functional lifespan.

The all-in-ONE Detergents deliver an average concentration that is 400% higher than the other enzyme detergents. One gallon can deliver 512 gallons (1938 liters) of cleaning and conditioning power. The dilution rate is .25 to .5 ounces per gallon (1.9 ml liters to 3.8 ml per liter) depending on the level of bioburden, encrustation, and stains.



all-in-ONE Detergent Enzyme Lubricant Cleaners,



easy FOAM-it Detergent Enzyme Lubricant Cleaners,

The all-in-ONE Enzyme Detergent Cleaner Lubricants are packaged as:

- UPS-32 easy-FOAM-it **Pre wash** Surgical Instrument Detergent Enzyme Lubricant Cleaners, 12 ea. 32 oz. per case,
- UPS-1 all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners, 4 gallons and 1 pump per case,
- UPS-5 all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners, 5 gallon jug,
- UPS-15 all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners, 15 gallon jug,
- UPS-30 all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners, 30 gallon jug,
- UPS-55 all-in-ONE Surgical Instrument Detergent Enzyme Lubricant Cleaners, 55 gallon jug,



Customer Service: 1-509-981-1553, email: Info@yourCEBA.com
Technical Support: 1-509-747-5027, email: Temple@yourCEBA.com