

your-CEBA Customer Service: 1-509-684-2457, email: RalphCEBA@live.com
Technical Support: John Temple, 1-509-747-5027, email: TempleCEBA@live.com

Endoscope Cleaning Sponges

Surgical Instrument Cleaning Sponges

Sponge Cleaning Kit for Medical Devices and Equipments

Sponges are available Dry and pre loaded with the all-in-ONE Enzyme Detergent

The Sponge is: soy based offering biodegradable benefits, retains more cleaning fluid than other sponges and paper based wipes, manufactured using RO purified diluents, and treated with UV Radiation prior to packaging. All Sponge products are periodically reviewed for integrity using laboratory replications and actual field storage conditions.

The all-in-ONE 4 Enzyme Detergent Cleaner is 400% times more effective than other enzyme cleaning solutions, delivering the four Enzymes (Protease, Amylase, Carbohydrase, and Lipase Enzymes) necessary to break down all forms of bioburden.

The all-in-ONE delivers the surface cleaning Surfactant Detergents needed for rapid and effective cleaning.

The all-in-ONE Enzyme Detergent Cleaners are non-irritating, biodegradable, neutral pH, and phosphate free.

Surgical Instrument Cleaning Sponges

D-Zyme Sponge, Order Number: ESW105200

Flat, Multi-Use, Sponge Loaded with

highly concentrated all-in-ONE Enzyme Detergents.

Individual Packaged - 200 ea. Per Case (4 boxes of 50)

Contents: Sponge Pre-soaked with Enzyme Detergent

Sponge Size: 4" Wide x 5" Long x 1/2" High

Order Number: ESW105200



Endoscope Cleaning Sponges

D-Zyme Sponge, Order Number: ESW104100

Tubular, Multi-Use, Sponge, Loaded with

highly concentrated all-in-ONE Enzyme Detergents.

Individual Packaged - 100 ea. Per Case (4 boxes of 25)

Contents: Sponge Pre-soaked with Enzyme Detergent

Sponge Size: 4.5" Long x 2.4" Diameter

Order Number: ESW104100

Sponge Cleaning Kit for Medical Devices and Equipment

Appli-Kit, Order Number: APT215035

Pouch contains Sponge and 215 ml of the highly concentrated all-in-ONE 4 Enzyme Detergents.

Pouch has a wide base pouch for stability and easy fluid transfer.

Individual Packaged - 35 ea. Per Case

Contents: Pouch with Sponge and 215 ml of Enzyme Detergent

Sponge Size: 4.5" Long x 2.4" Diameter

Order Number: APT215035

